




PATIENT SAFETY IS NOT NEGOTIABLE

DualCap[®]
Disinfection & Protection System

**DTXPlus[®]
& Safedraw[®]**
Disposable Transducer and
Closed in-line Sampling System



REDUCING MEDICAL ERRORS

DualCap® and DTXPlus® with Safedraw® products – engineered with several innovative features designed to prevent medical errors.

Straight Fluid Column

Clear housing easily identifies air in the line which can lead to damping and inaccurate information that guides decision making.


Safedraw Sample Port

No dead space, minimizes air, damping and inaccurate information that can result in errors.

The DualCap System

Locks In Protection

The only cap for male luer connectors that completely blocks disinfectant from the patient's bloodstream¹.



Safedraw Syringe

Closed draw system to reduce line contamination and unwanted air, critical features designed to reduce medical errors.

EasyVent® Luer Lock

Zero and prime without touch contamination assisting in the effort to protect against infection.

Angled Cannula

Minimizes air, reduces damping designed for improved accuracy to assist clinicians taking action.

Catalog Number	Description	Length	Tubing Color	Units Per Box/Carton
SAFEDRAW® CLOSED IN-LINE SAMPLING SYSTEM				
682142	Safedraw® with Meritran® DTXPlus® 60" Pressure Monitoring Kit and Blood Sampling System	120 / 30 cm	Clear	5/25
682143	Safedraw® with Meritran® DTXPlus® 72" Pressure Monitoring Kit and Blood Sampling System	150 / 30 cm	Clear	5/25
SAFEDRAW® STAND-ALONE SETS				
682144	Safedraw® 60" Blood Sampling System	120/30 cm	Clear	5/50
682145	Safedraw® 72" Blood Sampling System	150/30 cm	Clear	5/50
SAFEDRAW® ACCESSORIES				
682158	Safedraw® for BD Vacutainer® Transfer Device	N/A	N/A	25/125
682159	Safedraw® for BD Luer-Lok™ Safe-Needle	N/A	N/A	50/500

Catalog Number	Product Description	Caps Per Box	Boxes Per Case	Total Caps Per Case
DUALCAP® DISINFECTION & PROTECTION SYSTEM				
001-300	DualCap Solo™ Cap with foil lid for needle-free valve, packed in dispenser box.	350	12	4200
001-400	DualCap Solo™ Cap with foil lid for male luer connector, packed in dispenser box.	200	12	2400
500-LB	DualCap® IV Pole Strips® For needle-free valves. 35 strips/10 caps per strip.	350	12	4200
600-DB	DualCap® IV Pole Strips® For male luer connectors. 35 strips/5 caps per strip.	175	6	1050

¹ Avomeen Analytical Services, Testing Report prepared for Catheter Connections, October 2013.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 91 1238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115