

A close-up, black and white photograph of an elephant's head, showing its eye, textured skin, and large tusks. The elephant is looking slightly to the left.

BUILT FOR STRENGTH

The logo for HERO GRAFT, featuring the word "HERO" in a stylized, bold font with a registered trademark symbol, and "GRAFT" in a smaller, bold font below it.

HERO[®]
GRAFT

Venous Outflow Component

- Integrated, crush resistant radiopaque marker band
- Tougher blend of silicone for improved pushability



HeRO Graft


HeRO Graft (Hemodialysis Reliable OutFlow) is the **ONLY** fully subcutaneous AV access solution clinically proven to maintain long-term access for hemodialysis patients with **central venous stenosis**.

Silicone-Coated Venous Outflow Component

- No venous anastomosis
- Reinforced 48 braid nitinol: kink & crush resistant
- Removable and replaceable
- Radiopaque band (at distal tip)

Treatment Algorithm

Failing AVF or AVG due to central venous stenosis



Catheter-dependent patients

Improved Features of the HERO1001VOC

- Tougher blend of silicone for improved pushability
- New and enhanced radiopaque marker band
 - Integrated
 - Crush resistant



Ordering Information

Product Code	Component	Diameter (ID)	Length
HERO1001VOC	Venous Outflow Component w/new marker band	5mm	40cm (customizable)
HERO1002	Arterial Graft Component	6mm (ePTFE); 6mm - 5mm (connector taper)	53cm (connector: 3cm)
HERO1003	Accessory Component Kit	N/A	N/A

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 91 1238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115