


# Meritrans® DTXPlus®

## Disposable Transducer Sets

### KEY FEATURES




#### GABARITH™ TESTED

- Scientific test to ensure waveform accuracy of the entire system.
- Meets standard of Association for the Advancement of Medical Instrumentation (AAMI)



#### R.O.S.E.®

- Eliminates resonance overshoot via consistent damping through a fixed resistance.
- Works on pressure tubing between 36"-48".



#### M STOPCOCK

- Selected sets include M stopcock for blood sampling.
- Fits any common blood sampling device.

## Meritrans® DTXPlus® Disposable Transducer Sets

PRODUCT CODE	NO. OF DISPOSABLE TRANSDUCERS	FLUSH/ FLOW RATE (ML/HR)	TOTAL PRESSURE LINE LENGTH	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	BOX/CASE CONTENT (UNITS)
--------------	-------------------------------	--------------------------	----------------------------	-------------------------	------------------------	--------------------------

### DISPOSABLE TRANSDUCER SETS FOR THE INVASIVE MONITORING OF HEMODYNAMIC PRESSURE

682002	1	3	150cm (120/30)	Clear	-	5/50
682009	1	3	175cm (150/25)	Red	With R.O.S.E®	5/20
682032	2	3	2 x 150cm (120/30)	Red/Blue	-	5/25
682037	2	3	2 x 150cm (120/30)	Red/Blue	With 3rd IV line for add-on set	5/25
682041	2	3	2 x 150cm (120/30)	Blue/Yellow	-	5/25
682046	3	3	3 x 150cm (120/30)	Red/Blue/Yellow	-	5/25
682048	2	3	2 x 150cm (120/30)	Red/Blue	Velcro strap for arm mount	5/20

### DTXPLUS® TRANSDUCER SET WITH M STOPCOCK

688989	1	3	180cm (150/30)	Clear	Single Admin Set. Single Transducer Set with Integrated 1 x M stopcock (Red)	5/50
688990	2	3	180cm (150/30) 180cm (150/30)	Red Blue	Single Admin Set. Double Transducer Set with Integrated 1 x M stopcock (Red)	5/25

## Pressure Accessories

PRODUCT CODE	MODEL	PRODUCT DESCRIPTION	BOX/CASE CONTENT (UNITS)
682275	TMM/TBG	TBG with Snap-On TMM Cable Organiser (P/Ns 682093, 682274)	5/50
682276	TMM/TBG2	TBG2 with Snap-On TMM Cable Organiser	5/50

More configurations are available, please contact your Merit local representatives for details.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and direction for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

Merit Medical Asia Co., Ltd.  
Room 901B, 9th Floor,  
Empire Centre, 68 Mody Road,  
Tsim Sha Tsui, Kowloon, Hong Kong  
+852 2207 6888

China  
+86 10 8561 0788

Hong Kong  
+852 2207 6888

Indonesia  
+62 21 2939 8794

Korea  
+82 2 2205 0707

Malaysia  
+60 3 5624 4413

Myanmar  
+65 6750 9207

Philippines  
+65 6750 9207

Singapore  
+65 6750 9207

Taiwan  
+886 2 2656 5125

Thailand  
+65 6750 9207

Vietnam  
+84 28 3520 2922

Merit.com