



ACCURACY
1 mL AT A TIME

Medallion IZOFF™
Incremental Syringe

ACCURATE MANUAL
MEDICATION DELIVERY



Medallion IZOFF™

Incremental Syringe

The Medallion IZOFF Syringe has ribs at each mL along the syringe barrel providing an audible and tactile indication that 1 mL of fluid has been injected. These ribs, at every 1 mL, help ensure accurate and precise medication delivery.

\$2 Billion

The increased hospital costs alone of preventable adverse drug events affecting inpatients are about \$2 billion¹

FACT

Audible click and tactile feel for instant feedback on correct injection

Extremely accurate delivery



1.5M

Patients in the U.S., including people in hospitals and long-term care settings are harmed by medication errors²

FACT

Catalog number	Barrel Size	Increment Volumes
MIZ10	10 mL	1 mL
MIZ10P	10 mL	1 mL

1. The costs of adverse drug events in hospitalized patients. Adverse Drug Events Prevention Study Group. Bates DW1, Spell N, Cullen DJ, Burdick E, Laird N, Petersen LA, Small SD, Sweitzer BJ, Leape LL. <https://www.ncbi.nlm.nih.gov/pubmed/9002493>

2. The Joint Commission. National Patient Safety Goals Effective January 1, 2015. http://www.jointcommission.org/assets/1/6/2015_NPSG_HAP.pdf

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

merit.com

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30
(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115