

EXPERIENCE
THE POWER
OF ROTATION™

basix **Tau**™
INFLATION DEVICE

LEADING. LISTENING. INNOVATING.™

Merit Medical continues to innovate while improving legacy products customers trust. It's just one of the reasons Merit is a leader in the inflation market.

The basixTAU provides the quality and reliability known to Merit inflation devices and includes a revolutionary design of an added fold-out handle reducing the rotational force by 67%* during inflations.



Standard hand placement.




Fold-out handle hand placement.

*Versus basixKompak; Data on file

Easy-to-read luminescent analog display for low-light procedures. A maximum pressure of 30 ATM gauge.

20mL device capable of producing a maximum pressure of 30 ATM.

Fold-out handle feature brings versatility, ease of use, and additional ergonomic benefits.



ORDERING INFORMATION

Catalog Number	Description	Quantity per Box
IN4T130	IN4T130 - Standalone	5
IN4T530	IN4T130, MAP152 (Access Plus™ + hemostasis valve, metal wire introducer, torque device)	5
IN4T952	IN4T130, MAP952 (FLO50™ hemostasis valve, metal wire introducer, torque device)	5
IN4TXR52	IN4T130, FLO40XR52 (FLO40XR™ hemostasis valve, metal wire introducer, torque device)	5

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 91 1238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115