



# DiamondTOUCH™

Syringe



High resolution digital gauge provides exceptional visibility and accuracy



Offers the convenience of one-handed preparation



Up to  
**500**  
PSI

Simplified bone cement delivery and balloon inflation.

# DiamondTOUCH™

## Syringe

### 500 PSI, 30 ML SYRINGE

The DiamondTOUCH combines the visibility and accuracy of a digital device with the convenience of an ergonomic one-handed prep handle.

### DIAMONDTOUCH FEATURES

- One button functionality for simplicity and ease of use
- Large dial face with oversized time and pressure digits to enhance visibility
- Clutch mechanism allows for immediate cement flow stoppage
- 36-inch delivery line
- Ergonomic one-handed prep handle

### SPEED AND PERFORMANCE

- **Quick release ergonomic handle** - Releases easily, even from high-pressure inflations or cement deliveries.
- **One-handed preparation** - Preparation made easy with unique one-handed air removal. Simply press the handle against a solid surface.

### ROBUST AND RELIABLE

Strong components provide expanded capability with the security and confidence in performance that comes from using Merit products.



### ORDERING INFORMATION

Part Number	Description
IN9VCF	DiamondTOUCH Syringe

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

[merit.com](http://merit.com)

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115