



DO MORE WITH LESS

Splittable sheaths that require less
break force and less insertion force.

Splittable Sheath
Introducers



HYDROPHILIC PRELUDE SNAP™ SPLITTABLE SHEATH INTRODUCER

The Hydrophilic Prelude SNAP offers the same advanced hub design as the original Prelude SNAP, requiring 66% less break force than its predecessor.* **A sleek hydrophilic sheath coating reduces insertion force by 50%.*** Smooth sheath-to-dilator transitions help to reduce vessel trauma and support smooth navigation during your procedures.

Improved Grip

Expanded wings with non-slip, ridged surfaces facilitate a tighter, more assured grip.

Robust Valve Design

Tight seal to protect against air embolism and blood loss.

Hydrophilic Coating

Covers 75% of distal sheath length for ease of insertion while maintaining stability.

Seamless Transition

Smooth sheath-to-dilator transition to reduce vessel trauma and promote ease of insertion.



Easier to Break

Requires 66% less break force than its predecessor.*

ORDERING INFORMATION

Sideport

PRODUCT NUMBER (KIT)	SIZE	SHEATH LENGTH	HUB COLOR
PLSH-1006	6F	13 cm	Green
PLSH-1007	7F	13 cm	Orange
PLSH-1008	8F	13 cm	Blue
PLSH-1008.5	8.5F	13 cm	Light Blue
PLSH-1009	9F	13 cm	Black
PLSH-1009.5	9.5F	13 cm	Gray
PLSH-1010	10F	13 cm	Fuchsia
PLSH-1010.5	10.5F	13 cm	Light Fuchsia
PLSH-1011	11F	13 cm	Yellow
PLSH-1012	12F	13 cm	Brown
PLSH-1012.5	12.5F	13 cm	Light Brown
PLSH-2506	6F	25 cm	Green
PLSH-2507	7F	25 cm	Orange
PLSH-2508	8F	25 cm	Blue
PLSH-2509	9F	25 cm	Black
PLSH-2510	10F	25 cm	Fuchsia
PLSH-2510.5	10.5F	25 cm	Light Fuchsia
PLSH-2511	11F	25 cm	Yellow

Packaged 5 units per box.

Hydrophilic Prelude Snap Kits

- Splittable Sheath
- Dilator
- J-tipped Guide Wire, 0.038" x 50 cm or 0.038" x 80 cm
- Merit Advance® Needle, 18G x 7 cm
- Syringe

PRELUDE SNAP™ SPLITTABLE SHEATH INTRODUCER

The Prelude SNAP™ is one of the most advanced peelable sheaths on the market. Its high-quality hub requires 66% less break force than its predecessor,* a feature designed to facilitate effortless removal, helping to ensure leads stay in place. Seamless sheath-to-dilator transitions aid in reducing vessel trauma and providing ease of insertion—supporting optimal outcomes for your procedures.

Improved Grip
Expanded wings with non-slip, ridged surfaces facilitate a tighter, more assured grip.

Robust Valve Design
Tight seal to protect against air embolism and blood loss.

Easier to Break
Requires 66% less break force than its predecessor.*

Seamless Transition
Smooth sheath-to-dilator transition to reduce vessel trauma and promote ease of insertion.

ORDERING INFORMATION

Sideport			
PRODUCT NUMBER (KIT)	SIZE	SHEATH LENGTH	HUB COLOR
PLS-1006	6F	13 cm	Green
PLS-1007	7F	13 cm	Orange
PLS-1008	8F	13 cm	Blue
PLS-1008.5	8.5F	13 cm	Light Blue
PLS-1009	9F	13 cm	Black
PLS-1009.5	9.5F	13 cm	Gray
PLS-1010	10F	13 cm	Fuchsia
PLS-1010.5	10.5F	13 cm	Light Fuchsia
PLS-1011	11F	13 cm	Yellow
PLS-1012	12F	13 cm	Brown
PLS-1012.5	12.5F	13 cm	Light Brown
PLS-2506	6F	25 cm	Green
PLS-2507	7F	25 cm	Orange
PLS-2508	8F	25 cm	Blue
PLS-2509	9F	25 cm	Black
PLS-2510	10F	25 cm	Fuchsia
PLS-2510.5	10.5F	25 cm	Light Fuchsia
PLS-2511	11F	25 cm	Yellow

Non-Sideport			
PRODUCT NUMBER (KIT)	SIZE	SHEATH LENGTH	HUB COLOR
PLSX-1006	6F	13 cm	Green
PLSX-1007	7F	13 cm	Orange
PLSX-1008	8F	13 cm	Blue
PLSX-1008.5	8.5F	13 cm	Light Blue
PLSX-1009	9F	13 cm	Black
PLSX-1009.5	9.5F	13 cm	Gray
PLSX-1010	10F	13 cm	Fuchsia
PLSX-1010.5	10.5F	13 cm	Light Fuchsia
PLSX-1012	12F	13 cm	Brown
PLSX-1012.5	12.5F	13 cm	Light Brown
PLSX-2506	6F	25 cm	Green
PLSX-2507	7F	25 cm	Orange
PLSX-2508	8F	25 cm	Blue
PLSX-2509	9F	25 cm	Black
PLSX-2510	10F	25 cm	Fuchsia


Prelude Snap Kits

- Splittable Sheath
- Dilator
- J-tipped Guide Wire, 0.038" x 50 cm or 0.038" x 80 cm
- Merit Advance® Needle, 18G x 7 cm
- Syringe

Packaged 5 units per box.

PRELUDE PRESTIGE™ SPLITTABLE SHEATH INTRODUCER

Enjoy the same ease of insertion and break force of Prelude SNAP technology, now with the Prelude Prestige. A non-valved splittable sheath, the Prestige is designed for the placement of pacing/defibrillator leads and catheters. An improved grip delivers exceptional sheath peeling, and a wide range of sizes are available to accommodate your diverse clinical needs.



Prelude Prestige Kits

- Splittable Sheath
- Dilator
- J-tipped Guide Wire, 0.038" x 50 cm or 0.038" x 80 cm
- Merit Advance® Needle, 18G x 7 cm
- Syringe

ORDERING INFORMATION

PRODUCT NUMBER (KIT)	PRODUCT NUMBER (SHEATH/ DILATOR)	SIZE	SHEATH LENGTH	HUB COLOR
PLP-1005	PLPS-1005	5F	13 cm	Gray
PLP-1006	PLPS-1006	6F	13 cm	Green
PLP-1007	PLPS-1007	7F	13 cm	Orange
PLP-1008	PLPS-1008	8F	13 cm	Blue
PLP-1008.5	PLPS-1008.5	8.5F	13 cm	Light Blue
PLP-1009	PLPS-1009	9F	13 cm	Black
PLP-1009.5	PLPS-1009.5	9.5F	13 cm	Gray
PLP-1010	PLPS-1010	10F	13 cm	Fuchsia
PLP-1010.5	PLPS-1010.5	10.5F	13 cm	Light Fuchsia
PLP-1011	PLPS-1011	11F	13 cm	Yellow
PLP-1012	PLPS-1012	12F	13 cm	Brown
PLP-1012.5	PLPS-1012.5	12.5F	13 cm	Light Brown
PLP-1013	PLPS-1013	13F	13 cm	Purple
PLP-1014	PLPS-1014	14F	13 cm	Red
PLP-1015	PLPS-1015	15F	13 cm	Gray
PLP-1016	PLPS-1016	16F	13 cm	Green
PLP-2505	PLPS-2505	5F	25 cm	Gray
PLP-2506	PLPS-2506	6F	25 cm	Green
PLP-2507	PLPS-2507	7F	25 cm	Orange
PLP-2508	PLPS-2508	8F	25 cm	Blue
PLP-2509	PLPS-2509	9F	25 cm	Black
PLP-2510	PLPS-2510	10F	25 cm	Fuchsia
PLP-2511	PLPS-2511	11F	25 cm	Yellow
PLP-2512	PLPS-2512	12F	25 cm	Brown
PLP-2513	PLPS-2513	13F	25 cm	Purple
PLP-2514	PLPS-2514	14F	25 cm	Red
PLP-2515	PLPS-2515	15F	25 cm	Gray
PLP-2516	PLPS-2516	16F	25 cm	Green

Packaged 5 units per box. | Comes in full kits or individual sets

*Data on file

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733