


ELATION 5™




NEW 5-Stage Dilation Balloon

*Decrease frequency of needing
2 balloons*

Additional size offering

- 5 mm diameter balloon
- 21 mm diameter balloon



ORDERING INFORMATION

FIXED WIRE DILATION BALLOON

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)
EX6	5-6-7-8-9	8	2.8	180
EX8	7-8-9-10-11	8	2.8	180
EX10	9-10-11-12-13	8	2.8	180
EX12	11-12-13.5-15-16	8	2.8	180
EX15	14-15-16.5-18-19	8	2.8	180
EX18	17-18-19-20- 21	8	2.8	180

WIREGUIDED DILATION BALLOON

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)
EPX6	5-6-7-8-9	5.5	2.8	180
EPX8	7-8-9-10-11	5.5	2.8	180
EPX10	9-10-11-12-13	5.5	2.8	180
EPX12	11-12-13.5-15-16	5.5	2.8	180
EPX15	14-15-16.5-18-19	5.5	2.8	180
EPX18	17-18-19-20- 21	5.5	2.8	180
EPCX6	5-6-7-8-9	5.5	2.8	240
EPCX8	7-8-9-10-11	5.5	2.8	240
EPCX10	9-10-11-12-13	5.5	2.8	240
EPCX12	11-12-13.5-15-16	5.5	2.8	240
EPCX15	14-15-16.5-18-19	5.5	2.8	240
EPCX18	17-18-19-20- 21	5.5	2.8	240

UNIQUE SIZES are noted with **red/bold** part numbers

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733