



GASTROENTEROLOGY



EXPANDING THE POSSIBILITIES™

Stricture Management

EndoMAXX®

MAXXWIRE®

ALIMAXX-ES™

Elation5™

Elation®

BIG60®

BIG60ALPHA™

basixCOMPAK™

Monarch®

Unique Clinical Solutions

VacLok®

Stricture Management

EndoMAXX[®]

FULLY COVERED ESOPHAGEAL STENT

Accurate one-handed deployment

Enhanced anti-migration struts


Virtually no foreshortening or elongation

Increased flexibility & luminal patency


Metal proximal & distal repositioning sutures



Anti-migration struts



PROVEN LOW MIGRATION*



Catalog Number	Stent Size Options		Delivery Device		Packaged
	Stent Diameter (mm)	Stent Length (mm)	Working Diameter (F)	Working Length (cm)	
MAXX-197	19	70	24	90	1 unit
MAXX-1910	19	100	24	90	1 unit
MAXX-1912	19	120	24	90	1 unit
MAXX-1915	19	150	24	90	1 unit
MAXX-237	23	70	24	90	1 unit
MAXX-2310	23	100	24	90	1 unit
MAXX-2312	23	120	24	90	1 unit
MAXX-2315	23	150	24	90	1 unit

* Dua et al. GIE. Efficacy and safety of a new fully covered self-expandable non-foreshortening metal esophageal stent. 2014


MAXXWIRE[®] GUIDEWIRE

Super stiff main body delivers support and stability

5 cm straight atraumatic tip

1 cm distal measurement and progress tracking marks

Catalog Number	O.D. (in)	Length (cm)	Packaged
80129-404	0.035	180	1 unit
80129-406	0.035	260	1 unit



Stricture Management

ALIMAXX-ESTM

FULLY COVERED ESOPHAGEAL STENT


Accurate one-handed deployment

Unique small diameter stent sizes

Enhanced anti-migration struts

Virtually no foreshortening or elongation

Flexible tapered tip



EXPANDING TREATMENT OPTIONS



**UNIQUE
SIZES**

Catalog Number	Stent Size Options		Delivery Device		Packaged
	Stent Diameter (mm)	Stent Length (mm)	Working Diameter (F)	Working Length (cm)	
80129-207	12	70	22	67	1 unit
80129-208	12	100	22	64	1 unit
80129-209	12	120	22	62	1 unit
80129-210	14	70	22	67	1 unit
80129-211	14	100	22	64	1 unit
80129-212	14	120	22	62	1 unit
80129-213	16	70	22	67	1 unit
80129-214	16	100	22	64	1 unit
80129-215	16	120	22	62	1 unit
80129-201	18	70	22	67	1 unit
80129-202	18	100	22	64	1 unit
80129-203	18	120	22	62	1 unit
80129-204	22	70	22	67	1 unit
80129-205	22	100	22	64	1 unit
80129-206	22	120	22	62	1 unit

Stricture Management

ELATION⁵

BALLOON DILATOR

Worlds first 5-stage gastrointestinal
dilation balloon

Unique 5 mm & 21 mm diameter balloons

Decreases need for multiple balloons

Rapid and complete deflation

Optimal balloon clarity

Accurate, reliable balloon sizing



DO MORE WITH LESS

FIXED WIRE DILATION BALLOON 5-STAGE

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
EX6	5-6-7-8-9	8	2.8	180	1 unit
EX8	7-8-9-10-11	8	2.8	180	1 unit
EX10	9-10-11-12-13	8	2.8	180	1 unit
EX12	11-12-13.5-15-16	8	2.8	180	1 unit
EX15	14-15-16.5-18-19	8	2.8	180	1 unit
EX18	17-18-19-20- 21	8	2.8	180	1 unit


UNIQUE SIZES are noted with **red/bold** part numbers

WIREGUIDED DILATION BALLOON 5-STAGE


Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
EPX6	5-6-7-8-9	5.5	2.8	180	1 unit
EPX8	7-8-9-10-11	5.5	2.8	180	1 unit
EPX10	9-10-11-12-13	5.5	2.8	180	1 unit
EPX12	11-12-13.5-15-16	5.5	2.8	180	1 unit
EPX15	14-15-16.5-18-19	5.5	2.8	180	1 unit
EPX18	17-18-19-20- 21	5.5	2.8	180	1 unit
EPCX6	5-6-7-8-9	5.5	2.8	240	1 unit
EPCX8	7-8-9-10-11	5.5	2.8	240	1 unit
EPCX10	9-10-11-12-13	5.5	2.8	240	1 unit
EPCX12	11-12-13.5-15-16	5.5	2.8	240	1 unit
EPCX15	14-15-16.5-18-19	5.5	2.8	240	1 unit
EPCX18	17-18-19-20- 21	5.5	2.8	240	1 unit

UNIQUE SIZES are noted with **red/bold** part numbers

RELATION ORDERING GUIDE



Simplify inventory management with one of our 3 SKU options.



Widest range of diameters on a single balloon

Stricture Management

ELATION[®]

BALLOON DILATOR

3 distinct diameters at 3 corresponding pressures

Accurate, reliable balloon sizing

Optimal balloon clarity

Robust catheter design

Rapid and complete deflation



EFFICIENCY IN DILATION

FIXED WIRE DILATION BALLOON 3-STAGE

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
E6	6-7-8	8	2.8	180	1 unit
E8	8-9-10	8	2.8	180	1 unit
E10	10-11-12	8	2.8	180	1 unit
E12	12-13.5-15	8	2.8	180	1 unit
E15	15-16.5-18	8	2.8	180	1 unit
E18	18-19-20	8	2.8	180	1 unit

WIREGUIDED DILATION BALLOON 3-STAGE

Catalog Number	Inflated O.D. (mm)	Balloon Length (cm)	Min. Working Channel (mm)	Working Length (cm)	Packaged
EPB6	6-7-8	5.5	2.8	180	1 unit
EPB8	8-9-10	5.5	2.8	180	1 unit
EPB10	10-11-12	5.5	2.8	180	1 unit
EPB12	12-13.5-15	5.5	2.8	180	1 unit
EPB15	15-16.5-18	5.5	2.8	180	1 unit
EB18	18-19-20	5.5	2.8	180	1 unit
EPCB6	6-7-8	5.5	2.8	240	1 unit
EPCB8	8-9-10	5.5	2.8	240	1 unit
EPCB10	10-11-12	5.5	2.8	240	1 unit
EPCB12	12-13.5-15	5.5	2.8	240	1 unit
EPCB15	15-16.5-18	5.5	2.8	240	1 unit
ECB18	18-19-20	5.5	2.8	240	1 unit

Stricture Management

BIG60[®]

ANALOG INFLATION DEVICE

- Convenient all-in-one design
- Accurate dilation and rapid deflation
- Completely disposable



Catalog Number	Pressure Gauge	Pressure Rating	Size (mL)	Packaged
ENDO-AN6012	Analog	12 atm	60	5 units

BIG60[™] ALPHA

ANALOG INFLATION DEVICE

- Improved Ergonomics
- Easier Pressure Relief
- Rubberized Grip
- One Handed Prep
- Unique All-In-One Design
- Angled Pressure Gauge
- 3-way Stopcock Included
- Completely Disposable



Catalog Number	Pressure Gauge	Pressure Rating	Size (mL)	Packaged
EN-AL6012	Analog	12 atm	60	5 units



basix COMPACT™

ANALOG BILIARY INFLATION DEVICE

.....
Analog pressure gauge

.....
Gold standard for high pressure dilation

Catalog Number	Pressure Gauge	Pressure Rating (atm)	Size (mL)	Packaged
ENDO-AN2030	Analog	30	20	5 units



MONARCH®

DIGITAL BILIARY INFLATION DEVICE

.....
Digital pressure gauge

.....
Lighted screen

.....
Self-starting timer

Catalog Number	Pressure Gauge	Pressure Rating (atm)	Size (mL)	Packaged
ENDO-DI2030	Digital	30	20	5 units

Unique Clinical Solutions

VACLOK®

VACUUM PRESSURE SYRINGE

Polycarbonate barrel creates a “smooth as glass” feel

Ensures visual clarity of its contents


Plunger design holds negative pressure

Multiple negative pressure settings



Catalog Number	Size (mL)	Pressure Settings	Packaged
ENDO-VAC20	20	4 positions	20/box
ENDO-VAC20-SC*	20	4 positions	20/box
ENDO-VAC30	30	6 positions	25/box
ENDO-VAC30-SC*	30	6 positions	25/box
ENDO-VAC60	60	6 positions	25/box

*includes 1-way stopcock





Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com/endoscopy

EXPANDING THE POSSIBILITIES™

Please consult the product labels and inserts for any indications, contraindications, potential complications, warnings, precautions and directions for use.
All trademarks are the property of their respective owners.