


A large, natural rock archway made of reddish-brown sandstone. The arch frames a view of a blue sky with white clouds. A small person is visible standing on a rock ledge to the right, providing a sense of scale. The foreground shows a rocky, uneven ground surface.

ACCESS YIELDS RESULTS

Mini Access Kits

MINI ACCESS KITS

Merit is proud to offer a wide range of mini access devices that support your peripheral vascular access needs throughout the body.



MERIT MAK™ | MINI ACCESS KITS

The Merit MAK™ facilitates placement of a 0.035" (0.89 mm) or 0.038" (0.97 mm) diameter guide wire into the vascular system. The MAK kits allow you to start small, minimizing bleeding and trauma to the surrounding tissue. Kits are available in 4F and 5F and include the coaxial introducer/dilator pair, 21G needle, and 0.018" (0.46 mm) guide wire. For your convenience, all kit components are available separately. Merit MAK guide wires are available in many configurations to meet specific clinical needs. Guide wires are packaged in a fully enclosed hoop to ensure they are protected and ready for use.

KIT INCLUDES:

- Coaxial introducer/dilator pair
- 21G needle
- Guide wire



Catalog Number	French Size	Length	Guide Wire Diameter x Length	Guide Wire Core/Tip	Needle Gauge x Length	Qty/Box
STANDARD 7 CM NEEDLE KITS						
MAK401	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK501	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK401N	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK501N	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK501N-15	5	15 cm	0.018" (0.46 mm) x 60 cm	NT/PT	21G x 7 cm	10
MAK405	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
MAK505	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
ECHO-ENHANCED 7 CM NEEDLE KITS						
MAK401-E	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK501-E	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK401N-E	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK501N-E	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK405-E	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
MAK505-E	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
ECHO-ENHANCED 4 CM NEEDLE KITS						
MAK401-PED	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 4 cm	10
MAK501-PED	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 4 cm	10
MAK401N-PED	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	10
MAK405-PED	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 4 cm	10
MAK501N-PED	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	10
MAK505-PED	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 4 cm	10

SS = Stainless Steel PT = Platinum Tip NT = Nitinol

MERIT S-MAK™ | STIFFENED MINI ACCESS KITS

The Merit S-MAK™ facilitates placement of a 0.035" (0.89 mm) or 0.038" (0.97 mm) diameter guide wire into the vascular system. The S-MAK kit allows you to start small, minimizing bleeding and trauma to the surrounding tissue and is especially useful for access through difficult scar tissue. Kits are available in 4F and 5F sizes. For your convenience, all kit components are also available in separately sterilized packaging.

KIT INCLUDES:

- Coaxial introducer/stiffened dilator pair
- 21G echo-enhanced needle
- Guide wire




Catalog Number	French Size	Length	Guide Wire Diameter x Length	Guide Wire Core/Tip	Needle Gauge x Length	Qty/Box
STIFFENED DILATOR KITS WITH ECHO-ENHANCED 7 CM NEEDLE						
S-MAK401	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
S-MAK501	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
S-MAK401N	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
S-MAK501N	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
S-MAK501N-15	5	15 cm	0.018" (0.46 mm) x 60 cm	NT/PT	21G x 7 cm	10
S-MAK405	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
S-MAK505	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
S-MAK501N15BT	5	15 cm	0.018" (0.46 mm) x 60 cm	NT/PT	21G x 7 cm	10

SS = Stainless Steel PT = Platinum Tip NT = Nitinol BT = Bow Tie Guide Wire Insertion Device

MERIT PAK™ | PEDAL ACCESS KITS

The Merit PAK™ includes everything needed to access arteries below the knee. PAK kit micro-needles are engineered to access the small pedal arteries while introducers and dilators provide smooth transitions. Guide wires are packaged in a fully enclosed hoop for protection.



MINI PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Introducer & dilator
- Nitinol/Platinum guide wire
- Hemostasis Valve Adapter

STIFFENED PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Introducer & stiffened dilator
- Nitinol/Platinum guide wire
- Hemostasis Valve Adapter

HYDROPHILIC PEDAL ACCESS KIT

- Echo-enhanced, clear hub needle
- Hydrophilic sheath introducer with hemostasis valve and a side arm
- Nitinol/Platinum guide wire

Number	Type	Introducer OD FR	Introducer ID FR	Length	Guide Wire Diameter x Length	Guide Wire Core/Tip	Needle Gauge x Length	Qty/Box
ECHO-ENHANCED 4 CM NEEDLE KITS								
PAKM4NPD4	Mini	4	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	5
PAKM5NPD4	Mini	5	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	5
PAKS4NPD4	Stiffened	4	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	5
PAKS5NPD4	Stiffened	5	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	5
PAKPHR47NPD4	Hydrophilic	6	4	7 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	5
ECHO-ENHANCED 7 CM NEEDLE KITS								
PAKM4NPD7	Mini	4	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	5
PAKM5NPD7	Mini	5	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	5
PAKS4NPD7	Stiffened	4	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	5
PAKS5NPD7	Stiffened	5	3	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	5
PAKPHR47NPD7	Hydrophilic	6	4	7 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	5

MERIT MAK™ & S-MAK™ ACCESSORIES



Catalog Number	French Size	Introducer Length	Guide Wire Diameter x Length	Guide Wire Core/Tip	Needle Gauge x Length	Qty/Box
----------------	-------------	-------------------	------------------------------	---------------------	-----------------------	---------

COAXIAL INTRODUCER/DILATOR PAIRS (INDIVIDUALLY POUCHED)

MAK4F-10	4	10 cm	-	-	-	5
MAK5F-10	5	10 cm	-	-	-	5
MAK5F-15	5	15 cm	-	-	-	5
S-MAK4F-10	4	10 cm	(with stiffened dilator)	-	-	5
S-MAK5F-10	5	10 cm	(with stiffened dilator)	-	-	5
S-MAK5F-15	5	15 cm	(with stiffened dilator)	-	-	5

GUIDE WIRES (CONTAINED IN A HOOP AND INDIVIDUALLY POUCHED)

MAK001	-	-	0.018" (0.46 mm) x 40 cm	SS/SS	-	10
MAK001N	-	-	0.018" (0.46 mm) x 40 cm	NT/PT	-	10
MAK001N-60	-	-	0.018" (0.46 mm) x 60 cm	NT/PT	-	10
MAK001N-80	-	-	0.018" (0.46 mm) x 80 cm	NT/PT	-	10
MAK005	-	-	0.018" (0.46 mm) x 40 cm	SS/PT	-	10
MAK005-60	-	-	0.018" (0.46 mm) x 60 cm	SS/PT	-	10

NEEDLES (INDIVIDUALLY POUCHED)

AD21T21W	-	-	-	-	21G x 2.5 cm	25
AD21T41W	-	-	-	-	21G x 4 cm	25
AD21T71W	-	-	-	-	21G x 7 cm	25
AD21T91W	-	-	-	-	21G x 9 cm	25
AD21T41W-E*	-	-	-	-	21G x 4 cm	25
AD21T71W-E*	-	-	-	-	21G x 7 cm	25
AD21T91W-E*	-	-	-	-	21G x 9 cm	25

SS = Stainless Steel PT = Platinum Tip NT = Nitinol *Echo-enhanced needle

Catalog Number	Maximum Guide Wire Size	Qty/Box
BT01	≤ 0.018"/0.46mm	20

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733