




SMOOTH IS FAST

basix Sky™
Inflation Device



Smooth Is Fast

The basixSKY™ Inflation Device is a +30 ATM/BAR (zero to +441 psi) advanced analog inflation device designed for fast inflation and easy one-handed preparation, streamlining your procedures when speed and precision are critical.



COMPLETE PROCEDURAL SOLUTIONS

The basixSKY Inflation Device is part of our commitment to offering complete solutions that address unique clinical needs.

The basixSKY with the PhD™
An innovative hemostasis valve with push-and-release bleedback seal technology.

ORDERING INFORMATION

Catalog Number	Product Description	Units Per Box
SKY1130	basixSKY™ Inflation Device	10
SKY1152	basixSKY™ Inflation Device with MAP152 – AccessPLUS™ Hemostasis Valve, metal insertion tool, torque device	10
SKY1302	basixSKY™ Inflation Device with MAP302 – Honor® Hemostasis Valve, metal insertion tool, torque device	10
SKY1802	basixSKY™ Inflation Device with MAP802 – PhD™ Hemostasis Valve, metal insertion tool, torque device	10

Intended Purpose/Use: The basixSKY Inflation Device is used to inflate and deflate an angioplasty balloon or other interventional device and to measure the pressure within the balloon.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733