

# AEROmini<sup>®</sup>

NEW  
SIZES!



 **MERTMEDICAL**  
ENDOTEK<sup>®</sup>

# AEROmini<sup>®</sup>

- Fully Covered
- Enhanced proprietary anti-migration struts
- Virtually no foreshortening or elongation
- Hydrophilic coating



Catalog Number	Stent Size Options		Delivery Device		Package
	Stent Diameter (mm)	Stent Length (mm)	Working Diameter (Fr)	Working Length (cm)	

## THROUGH-THE-SCOPE

<b>NEW</b>	<b>AERO610</b>	6	10	7.9	110	1 unit
<b>NEW</b>	<b>AERO615</b>	6	15	7.9	110	1 unit

## OVER-THE-WIRE

<b>NEW</b>	<b>AERO810</b>	8	10	12.5	70	1 unit
	<b>AERO815</b>	8	15	12.5	70	1 unit
	<b>AERO820</b>	8	20	12.5	70	1 unit
<b>NEW</b>	<b>AERO1010</b>	10	10	12.5	70	1 unit
	<b>AERO1015</b>	10	15	12.5	70	1 unit
<b>NEW</b>	<b>AERO1210</b>	12	10	16	70	1 unit
	<b>AERO1215</b>	12	15	16	70	1 unit
<b>NEW</b>	<b>AERO1410</b>	14	10	16	70	1 unit
	<b>AERO1415</b>	14	15	16	70	1 unit

Indications for Use: The MERIT ENDOTEK AEROmini Tracheobronchial Stent System is indicated for use in the treatment of tracheobronchial strictures produced by malignant neoplasms.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733