



DISCOVER DISTAL

PreludeSYNC^oDISTAL™

Radial Compression Device



THE DISTAL DIFFERENCE

Distal access brings versatility to interventional procedures as an additional radial access site option that can offer benefits for both the patient and operator.

Secure

Three-strap band stays securely in place and is specially designed for distal radial artery access.

Comfortable

Unique band configuration is soft, comfortable, and adjustable for larger hands.

Effective

Provides clear, unobstructed site visibility and focused compression.

DISTAL BENEFITS:

- Shorter hemostasis time*
- Lower incidence of radial artery occlusion*
- Operator ergonomics
- Patient comfort
- An additional radial access site option

Optional extension band (included)



Slip-n-Lock™ cap

- Specialized connection¹ for safety during inflation.
- Secures syringe during inflation/deflation.
- If a PreludeSYNC DISTAL syringe is not available, remove the Slip-n-Lock cap² to allow deflation with a standard luer syringe.

ORDERING INFORMATION

Distal Compression Device with Syringe	PreludeSYNC Distal Compression Device Description	Units/Box
SDRB-REG-RT	Right Hand Device (Specialized Connection)	5
SDRB-REG-LT	Left Hand Device (Specialized Connection)	5

* Kiemeneij F et al. 2019. "Keep an Open Mind about Distal Radial Access." *Cath Lab Digest* 27, no. 3 (March 5). <https://www.cathlabdigest.com/content/keep-open-mind-about-distal-radial-access>

1. Patent pending. 2. CAUTION: DO NOT remove the Slip-n-Lock connector while sheath is still in the patient.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733