



POWER & CONTROL



CCS[™]

CORONARY CONTROL SYRINGES

IT ALL BEGAN WITH THE CORONARY CONTROL SYRINGE

Since 1987 Merit Medical has been manufacturing and selling a variety of coronary control syringes. With more than 30 years of experience, Merit has designed and created quality control syringes with features that continue to meet ever-changing clinical expectations.



Clear Barrel

Provides exceptional clarity and smooth-as-glass feel

Solid Plunger Body

Maintains stability and durability under pressure



Rotating or Fixed Luers

Allow for flexibility and confidence in luer connection

Safety Space

Intended to minimize the potential for air bubbles introduced into the catheter



Smart Tip Option

Creates a unique and smooth feel during injection

Plunger/Grip Options

Multiple configurations for plungers and grips on every syringe



Visit **Merit.com** or contact your Account Executive to find out how Merit can help you deliver efficient, patient-focused, and cost-effective care.

INJECT 8® Inject 10ⁿ®

The evolution of the coronary control syringe has now provided two revolutionary syringe designs. The Inject8™ and Inject10n™ allow clinicians to generate more pressure with less applied force.



Improved Visualization

Smaller ID barrels improve visualization with less contrast

Increased Flow

Smaller ID barrels increase flow through smaller french-sized catheters

Decreased Hand Fatigue

Smaller ID barrels require less force and reduce hand fatigue during injections

CORONARY CONTROL SYRINGE ORDERING INFORMATION

Catalog Number	Size (Volume)	Hub Style	Grip	Plunger Style
CCS100	12 mL	Rotating Male Luer Adapter	Ring Grip	Locking Ring Plunger
CCS101	12 mL	Fixed Male Luer Connector	Ring Grip	Locking Ring Plunger
CCS200	12 mL	Rotating Male Luer Adapter	Ring Grip	Solid Palm Pad Plunger
CCS201	12 mL	Fixed Male Luer Connector	Ring Grip	Solid Palm Pad Plunger
CCS320	20 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
CCS321	20 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger
CCS460	6 mL	Rotating Male Luer Adapter	Ring Grip	Locking Ring Plunger
CCS461	6 mL	Fixed Male Luer Connector	Ring Grip	Locking Ring Plunger
CCS600	12 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
CCS601	12 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger
CCS800	10 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
CCS801	10 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger
CCSB700	12 mL	Rotating Male Luer Adapter	Ring Grip	Solid Palm Pad Plunger
CCSB701	12 mL	Fixed Male Luer Connector	Ring Grip	Solid Palm Pad Plunger
CCSB710	12 mL	Rotating Male Luer Adapter	Wing Grip	Solid Palm Pad Plunger
CCSB711	12 mL	Fixed Male Luer Connector	Wing Grip	Solid Palm Pad Plunger
CCSB900	10 mL	Rotating Male Luer Adapter	Ring Grip	Solid Palm Pad Plunger
CCSB901	10 mL	Fixed Male Luer Connector	Ring Grip	Solid Palm Pad Plunger
CCSB910	10 mL	Rotating Male Luer Adapter	Wing Grip	Solid Palm Pad Plunger
CCSB911	10 mL	Fixed Male Luer Connector	Wing Grip	Solid Palm Pad Plunger

Smart Tip™

ST800	10 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
ST801	10 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger
ST600	12 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
ST601	12 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger
ST320	20 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
ST321	20 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger

Inject8™

CCS880	8 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
CCS881	8 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger
CCSB880	8 mL	Rotating Male Luer Adapter	Ring Grip	Solid Palm Pad Plunger
CCSB881	8 mL	Fixed Male Luer Connector	Ring Grip	Solid Palm Pad Plunger
CCSW880	8 mL	Rotating Male Luer Adapter	Wing Grip	Solid Palm Pad Plunger
CCSW881	8 mL	Fixed Male Luer Connector	Wing Grip	Solid Palm Pad Plunger

Inject10n™

CCX010	10 mL	Rotating Male Luer Adapter	Ring Grip	Solid Ring Plunger
CCX011	10 mL	Fixed Male Luer Connector	Ring Grip	Solid Ring Plunger
CCXB010	10 mL	Rotating Male Luer Adapter	Ring Grip	Solid Palm Pad Plunger
CCXB011	10 mL	Fixed Male Luer Connector	Ring Grip	Solid Palm Pad Plunger

Packaged 25 units per box, 4 boxes per case.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733