

ReSolve MINI™

Locking Drainage Catheter



5

10



**NEW
PRODUCT**

14F

20



merit.com

The ReResolve Mini™ Locking Drainage Catheter pigtail is up to 42% smaller than our standard locking pigtails, allowing the catheter to be placed into smaller fluid pockets. It features a tapered tip and hydrophilic coating for a smooth introduction, designed to reduce patient discomfort.

Multiple large holes
Help maximize drainage

Atraumatic
Internal pigtail designed to reduce tissue irritation

Radiopaque marker band
Denotes where the pigtail is to ensure all of the holes are in the fluid pocket

Depth markers
Located on the catheter shaft aid accurate placement and act as a reference for correct positioning

42%
SMALLER



ReResolve Locking



ReResolve Mini

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.

Merit Medical Systems, Inc. 1600 West Merit Parkway • South Jordan, Utah 84095

Main: +1.801.253.1600 Customer Service: +1.801.208.4300

406385001_001 ID120219