


# Cultura™ FLEX SWAB

Designed for patient comfort, the Cultura Flex Swab provides superior collection and release of specimens. Available in two breakpoint configurations.



**Soft Flocked Tip,  
16 mm Length**  
Short fibers create a thin absorbent layer for optimal sample collection.

**Flexible Shaft,  
150 mm Overall Length**  
Ergonomic design improves patient comfort.

**Swab Breakpoint**  
Indicator mark on handle for easy break off of swabstick into transport medium.

Catalog Number	Description	Packaged
C100F	100 mm breakpoint, 150 mm overall length	250 per box
C80F	80 mm breakpoint, 150 mm overall length	250 per box