

# Cultura™ TRANSPORT SYSTEM

The Cultura Transport System contains Viral Transport Medium specially formulated for virus preservation, offering the flexibility of either room temperature or refrigerated transport.

## VIRAL TRANSPORT MEDIUM, 3 mL

- **Temperature Range Flexibility** Room temperature or refrigerated (up to 120 hours) after specimen collection.
- **Maintains Sample Integrity** Formulated with proteins for virus preservation and antibiotics to inhibit bacteria and fungi overgrowth.
- **pH Specific** Phenol red pH indicator ensures medium integrity. Buffer solution maintains neutral pH.
- **Reduces Risk of Biohazard Spills** Self-standing, flat-bottom conical tube.
- **Prevents Leakage** Twist-on cap with dual-sealing surface design.
- **Traceable** Unique QR code on each vial enables efficient tracking of samples in Laboratory Information Management System software.



Catalog Number	Description	Packaged
VDCV100	1-Vial with 3 mL viral transport medium	100 per box

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.