



basix **ALPHA**TM Inflation Device



Built for Better Inflation

The basixALPHATM Inflation Device is an analog inflation syringe with a +30 ATM/BAR (zero to + 441 PSI) capability and 20 mL capacity, designed for use in a variety of interventional procedures. More than just an accessory, the basixALPHA has a number of product features which increase ease of use and solve common clinical challenges.



Unobstructed Visibility

A clear polycarbonate barrel allows for better visualization of bubbles in the device.

Fast Inflation Times, Minimal Effort

Coarse threads and a smaller diameter barrel enable quick inflation times. Rotational torque and number of revolutions to reach pressure have been minimized with the intent to reduce clinician hand fatigue.



The basixALPHA with the PhDTM

An innovative hemostasis valve with push-and-release bleedback seal technology.

Simplified One-Handed Prep Touch

Depress the handle for one-handed prep, freeing the opposite hand.

Ergonomic Handle

Engineered to fit the natural curves and sizes of operators' hands.

COMPLETE PROCEDURAL SOLUTIONS

The basixALPHA Inflation Device is part of Merit Medical's commitment to offering complete solutions that address unique clinical needs.

ORDERING INFORMATION

Catalog Number	Product Description	Units Per Box
IN3130	basixALPHA TM Inflation Device	5
IN3152	basixALPHA TM Inflation Device with MAP152 – AccessPLUS TM hemostasis valve, metal insertion tool, torque device	5
IN3302	basixALPHA TM Inflation Device with MAP302 – Honor [®] hemostasis valve, metal insertion tool, torque device	5
IN3802	basixALPHA TM Inflation Device with MAP802 – PhD TM hemostasis valve, metal insertion tool, torque device	5

Indications for Use: The basixALPHA Inflation Device is used to inflate and deflate an angioplasty balloon or other interventional device and to measure the pressure within the balloon.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.[™]

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733