

BIG 60 ALPHA™

Built for Better Inflation




**IMPROVED
ERGONOMICS**

**MERTMEDICAL
ENDOTEK®**

BIG60 ALPHA™

Built for better inflation

- Improved Ergonomics
- Easier Pressure Relief
- Rubberized Grip
- One Handed Prep
- Unique All-In-One Design
- Angled Pressure Gauge
- 3-way Stopcock Included
- Completely Disposable



Catalog Number	Pressure Gauge	Pressure Rating	Volume	Packaged
EN-AL6012	Analog	12 ATM	60 mL	5 units

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com/endoscopy

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
Main: +1.801.253.1600
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle
East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733