



SCOUT Bx™

Biopsy Delivery System

LOCALIZE AT THE TIME OF BIOPSY

INTEGRATED APPROACH TO LOCALIZATION WITH SCOUT

- Eliminates secondary localization procedures
- Helps reduce time to surgery
- Minimizes patient's trips to the facility
- End deploy designed to work with the most popular stereotactic and MRI guided biopsy devices
- Cleared for long-term implant with no restriction on the length of time the reflector can remain implanted

PRECISE DELIVERY SYSTEM

Compatible with:

- Hologic Eviva 0913-20, 1213-20 (STX)
- Hologic BREV09 (STX)
- ATEC ILS 0914-20 (MRI)



ORDERING INFORMATION

Model Number	Product Description
SSR13B-01	SCOUT Bx Reflector and Delivery System End Deploy 15G 12.3 cm

SIZE REFERENCE



To order, call 866.920.9444

Questions? Please contact your local Merit Oncology Representative or email our Customer Service Team at cs@merit.com

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
 1600 West Merit Parkway
 South Jordan, Utah 84095
 1.801.253.1600
 1.800.35.MERIT

Merit Medical Europe, Middle
 East & Africa (EMEA)
 Amerikalaan 42, 6199 AE
 Maastricht-Airport
 The Netherlands
 +31 43 358 82 22

Merit Medical Ireland Ltd.
 Parkmore Business Park West
 Galway, Ireland
 +353 (0) 91 703 733