

*Aspira*<sup>®</sup>  
DRAINAGE SYSTEM

# PHYSICIAN CHOICE OF BAG WITH SIPHON ASSIST OR BOTTLE



## FASTER MORE-COMFORTABLE DRAINAGE COMPARED TO PLEURX

- Saves patient time
- A faster drainage flow rate means patients can get back to what they enjoy
- 67% more catheter fenestrations<sup>1</sup>
- 72% larger valve Inner Diameter<sup>2</sup>
- Siphon assist on bag
- Less painful/discomfort
- The siphon-assist design of the Aspira Bag produces a consistent pressure flow for minimized patient discomfort

## PATIENT CONVENIENCE

- Ease of use for patient
- Aspira Bags are more compact for travel
- Requires less storage space at home

**68%** higher average flow rate<sup>1,3,4</sup>

**23%** less storage space<sup>1</sup>


**Table 1.** Simulated average flow rate comparison<sup>1,3,4</sup>

Aspira Drainage Bag	Aspira Drainage Bottle	PleurX Drainage System
7.46 mL/sec	4.76 mL/sec	4.45 mL/sec

Twenty Aspira Drainage Bottles and dressings take up 23% less storage space than 20 PleurX Bottles and dressings.

**71%** less storage space<sup>5</sup>

Twenty Aspira Drainage Bags and dressings require approximately the same amount of space as one box of four PleurX Bottles and dressings, reducing storage space and waste.



“There was just no way I was going down without a fight because I could have years ahead of me yet... It’s got to be under my conditions. Every time we get a delivery of **Aspira, it’s more life, more freedom, more adventures.**”

—Andrea Peterson,  
Aspira Patient

## FEEL THE DIFFERENCE IN CATHETER PLACEMENT



*Aspira Drainage System Tunneler*

### ENHANCED TUNNELER DESIGN

For placement flexibility, Aspira's malleable, pre-bent tunneler is 37% longer than the PleurX.<sup>TM1</sup>



*Prelude SNAP™ 165 Valved Splittable Sheath Introducer*

### VALVED INTRODUCER

A peel-away introducer sheath helps you reduce fluid leakage and/or air aspiration during placement.<sup>1</sup>



### REPLACEABLE VALVE

The valve can be replaced without the need for a catheter exchange.

## WHITE-GLOVE SERVICE



We offer a 24-hour Aspira Support Line for patients and healthcare professionals.

**1-833-3ASPIRA**  
(1-833-327-7472) or

**1-800-35MERIT**  
(1-800-356-3748)



Find physician and patient resources and supplies at [myAspira.com](http://myAspira.com).



Our support team provides hospital planning and product education to simplify and ease transitions.

## Aspira kit components



### INSERTION TRAY COMPONENTS

Catalog Number: Peritoneal 4992206, Pleural 4992506

- |   |  |
|---|--|
| 1 – Drainage catheter with stylet, 15.5F                        | 1 – Slide clamp                          |
| 1 – #11 Futura® Safety Scalpel                                  | 1 – Aspira valve assembly                |
| 1 – Introducer needle, 18G                                      | 1 – Aspira valve cap                     |
| 1 – J-tip guide wire, 0.038" x 80 cm                            | 1 – Suture wing                          |
| 1 – Dilator, 8F   | 1 – Luer adapter                         |
| 1 – Dilator, 12F  | 1 – Universal tubing adapter             |
| 1 – Prelude SNAP 165 Valved Splittable Sheath Introducer, 16.5F | 2 – Split gauze, 4" x 4" (10 cm x 10 cm) |
| 1 – Tunneler  | 6 – Gauze, 4" x 4" (10 cm x 10 cm)       |
|   | 1 – Surgical tape                        |
|   | 1 – Adhesive dressing                    |



### DRESSING KIT COMPONENTS

Catalog Number: 4991503

- |  |                         |
|--|-------------------------|
| 1 – Adhesive dressing                    | 1 – Tape strips (3 per) |
| 2 – Split gauze, 4" x 4" (10 cm x 10 cm) | 3 – Alcohol wipes       |
| 4 – Gauze, 4" x 4" (10 cm x 10 cm)       | 1 – CSR wrap            |
| 1 – Gloves (1 pair)                      |                         |



### DRAINAGE KIT COMPONENTS (BAG)

Catalog Number: 4992301

- |                             |                         |
|-----------------------------|-------------------------|
| 1 – Drainage bag (1,000 mL) | 1 – Tape strips (3 per) |
| 2 – Alcohol wipes           | 1 – Aspira valve cap    |
| 1 – Slide clamp             |                         |



### DRAINAGE KIT/DRESSING KIT COMPONENTS (BOTTLE)

Catalog Number: ASP-4 (4 PACK), ASP-10 (10 PACK)

- |                                |  |
|--------------------------------|--|
| 1 – Drainage bottle (1,000 mL) | 2 – Split gauze, 4" x 4" (10 cm x 10 cm) |
| 3 – Alcohol wipes              | 4 – Gauze, 4" x 4" (10 cm x 10 cm)       |
| 1 – Slide clamp                | 1 – Gloves (1 pair)                      |
| 1 – Tape strips (3 per)        | 1 – CSR wrap                             |
| 1 – Aspira valve cap           |  |
| 1 – Adhesive dressing          |  |



### LUER ADAPTER COMPONENTS

Catalog Number: 4992305

- 1 – Luer adapter
- 1 – Aspira valve cap



### UNIVERSAL TUBING ADAPTER COMPONENTS

Catalog Number: 4992307

- 1 – Universal tubing adapter
- 1 – Aspira valve cap



### VALVE REPAIR KIT COMPONENTS

Catalog Number: 4992306

- 1 – Slide clamp
- 1 – Aspira valve assembly
- 1 – Aspira valve cap

## Accessories kit components



### CATHETER INSERTION ACCESSORIES KIT

Catalog Number: K20T-00033

- 1 – Pre-filled syringe saline, 10 mL
- 1 – Medallion® Syringe, 10 mL
- 1 – Bag with adhesive, 6" x 10"
- 2 – ChloroPrep®, 10.5 mL (tinted)
- 1 – Suture silk 2-0, 30"
- 1 – Suture VICRYL® 4-0, 18"
- 1 – Lidocaine 1%, 50 mL
- 1 – ShortStop® Temporary Sharps Holder
- 1 – Vial access device, 20 mm
- 1 – Drape, 72" x 118" with 4" x 6" fenestration
- 1 – Alcohol prep pad
- 1 – Safety needle, 22G x 1 1/2"
- 1 – Safety needle, 25G x 1"
- 1 – Mosquito forceps straight, 4 3/4"
- 1 – Adson forceps, 4 3/4"
- 1 – Tray, 9" x 5" x 1"
- 1 – CSR wrap, 36" x 36"

## ORDERING INFORMATION

### Aspira Drainage System

Peritoneal (15 cm from cuff to proximal fenestration)	Pleural (11 cm from cuff to proximal fenestration)	Description
4992206	4992506	<b>Catheter Kit</b> 1 Insertion tray with catheter

### Aspira Drainage System with Bag

Peritoneal (15 cm from cuff to proximal fenestration)	Pleural (11 cm from cuff to proximal fenestration)	Description
4992207	4992507	<b>Complete Kit</b> 1 Insertion tray with catheter, 5 drainage kits, 5 dressing kits, voucher for 20 drainage kits and 20 dressing kits
4992208	4992508	<b>Essential Kit</b> 1 Insertion tray with catheter, 20 drainage kits, 20 dressing kits
4992209	4992509	<b>Starter Kit</b> 1 Insertion tray with catheter, 5 drainage kits, 5 dressing kits



Aspira  
Drainage  
System  
with bag

### Aspira Drainage System with Bottle

Peritoneal (15 cm from cuff to proximal fenestration)	Pleural (11 cm from cuff to proximal fenestration)	Description
4992204	4992504	<b>Starter Kit</b> 1 Insertion tray with catheter, 4 drainage kits/dressing kits
4992224	4992524	<b>Complete Kit</b> 1 Insertion tray with catheter, 4 drainage kits/dressing kits, voucher for 20 drainage kits/dressing kits



Aspira  
Drainage  
System  
with  
bottle

### Aspira Drainage Accessories

Catalog Number	Description
4991503*	Dressing kit
4992301*	Drainage kit with bag
ASP-4**	Drainage kit with bottle/dressing kit
ASP-10***	Drainage kit with bottle/dressing kit
4992305*	Luer adapter
4992306*	Valve repair kit
4992307*	Universal tubing adapter

\* Packaged 5 units per box \*\* Packaged 4 units per box \*\*\* Packaged 10 units per box

Visit [myaspira.com](http://myaspira.com) for more information.

All trademarks are the property of their respective owners.

1. Data on File.
2. Data on file. Internal diameters measured at the narrowest point. N=10 per device.
3. Data on file. Simulated water drainage. N=9 for Aspira Drainage Bag, N=8 for PleurX Drainage Bottle. Bench data may not represent clinical outcomes. Different test methods may yield different results.
4. Data on file. Simulated water drainage under similar conditions. N=40 for Aspira Drainage Bag, N=8 for PleurX Drainage Bottle. Bench data may not represent clinical outcomes. Different test methods may yield different results.
5. Data on file. PleurX packaging measurements represent five boxes of PleurX Starter Kits containing four 1,000-mL drainage bottles and dressings per box. Aspira packaging measurements represent one box of an Aspira Essential Kit containing twenty 1,000-mL drainage bags and dressings per box.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.

### Catheter Insertion Accessories Kit

Catalog Number	Description
K20T-00033	Contains additional accessories to help with catheter insertion

Packaged 6 units per box

### Prelude SNAP165 Valved Splittable Sheath Introducer

Catalog Number	Description
PLS165	Valved splittable sheath introducer, 16.5F

Packaged 10 units per box. Requires 6–8 weeks lead time.



Understand. Innovate. Deliver.™

[merit.com](http://merit.com) | [myaspira.com](http://myaspira.com)

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733