

# SAFEGUARD FOCUS™

## COMPRESSION DEVICE

Achieve targeted hemostasis with the SafeGuard Focus,™ a specialized compression device used to protect closed surgical sites.\*

Available in two unique applications—a contoured adhesive dressing and an adhesive-free flexible strap—the SafeGuard Focus is specifically designed for use in patients with pacemakers and implantable cardioverter defibrillator (ICD) pockets.

### ADHESIVE



### ADHESIVE-FREE



#### Cooling

Available with cooling in both adhesive and adhesive-free configurations.

#### Unobstructed Site Visibility

Clear window facilitates surgical site assessment without device removal.

#### Comfortable Design

Soft thermoplastic balloon and nylon loop enhance patient comfort.

#### Durable Attachment

Medical grade adhesive backing or VELCRO® strap offer secure attachment.\*\*

#### Ease of Use

Quick application without adjusting patient position.\*\*\*

Catalog Number	Description
<b>SGF1A</b>	Contoured adhesive dressing
<b>SGF1AC</b>	Contoured adhesive dressing with cooling
<b>SGF1V</b>	Adhesive-free flexible strap
<b>SGF1VC</b>	Adhesive-free flexible strap with cooling

Packaged 5 units per box

\*Not indicated for femoral artery compression. | \*\*The adhesive portion of the SafeGuard Focus device should not be used over damaged skin. | \*\*\*Adhesive version only. VELCRO® is a registered trademark of Velcro BVBA. Used with permission

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.

### SAFEGUARD FOCUS Cool™

#### COMPRESSION DEVICE

By adding saline to a pre-loaded syringe, the SafeGuard Focus Cool applies cooling directly to the application site while delivering compression.