



PreludeSYNC^oEZTM

Radial Compression Device

A close-up photograph showing a human hand being held by a tan-colored medical mannequin hand. A clear plastic device, the PreludeSYNC EZ, is being applied to the wrist. The device has a white strap and a green base with the 'MERIT EZ' logo. The device is positioned over the radial artery. The background is a light blue surface.

**CONTROLLED
PLACEMENT.
PRECISE
COMPRESSION.**

The clear choice for hemostasis.

Controlled Placement. Precise Compression.

Be confident in your radial interventions with the PreludeSYNC EZ, an advanced solution for achieving patent hemostasis. Clear visibility provides quick detection and precise compression while a unique cinch-style band offers controlled two-handed placement. The PreludeSYNC EZ simplifies application for you and your staff without sacrificing patient comfort.



Hemostasis in a Cinch

Special cinch-style band allows for excellent control during placement.

Clarity Is Confidence

- Large clear window supports accurate placement over arteriotomy and early detection of hematomas.
- Reduced fill volume facilitates quick deflation.

Safety & Convenience

Safeguard your inflations with the secure Slip-n-lock[™] cap while still preserving deflation options using a standardized syringe.



Inflate and deflate with patented Slip-n-Lock connection



Optional deflation with standard syringe*

Watch the PreludeSYNC EZ in action at merit.com/preludesynceez

ORDERING INFORMATION

Catalog Number	Product Description
EZ-SML	26 cm radial compression device
EZ-REG	32 cm radial compression device

Packaged 5 units per box. | Not made with natural rubber latex.

*CAUTION: DO NOT remove the Slip-n-Lock cap while sheath is still in the patient.

**Kink resistance versus leading competitor. 6F data on file.

Complete Radial Solution

Prelude IDEAL[™] • Braided thin-wall technology
Hydrophilic Sheath Introducer • 23% more kink resistant**

To learn more, speak with your Merit Sales Representative or visit merit.com/IDEAL.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733