



DuraFlow[™]

Pre-curve Chronic Hemodialysis Catheter

DuraFlow™

Pre-curve Chronic Hemodialysis Catheter

Engineered for Durability

15.5 F DuraFlow™ catheter

Facilitates high flow rates at modest arterial pressures. Lumen design provides improved kink resistance and flexibility in challenging tunnel locations.

Durathane® catheter material affords strength, softness, and resistance to a wide range of commonly used site care agents.




EmboSafe™

Peel-away valved sheath provides safety and security of vascular access and greatly reduces the risk of air embolism.



Oval guidewire holes allow for a smoother wire transition.

Robust polyester cuff material allows optimal tissue ingrowth and catheter securement.



Rotating suture wings ensure patient comfort and flexibility to position catheter.

ORDERING INFORMATION

Catalog Number	Catheter Length
BASIC KIT	
-	20 cm
-	22 cm
H787103012065	24 cm
H787103012075	28 cm
H787103012085	32 cm
H787103012095	36 cm

Packaged 5 units per box. | Catheter measurements from tip to hub.

BASIC KIT INCLUDES:

- 1-DuraFlow Hemodialysis Catheter
- 1-16F Single Valve Peelable Introducer Sheath
- 1-12F Vessel Dilator
- 1-14F Vessel Dilator
- 1-Tri-Ball Tunneler with Sleeve
- 2-Injection Caps
- 1-Adhesive Wound Dressing
- 1-#11 Safety Scalpel
- 1-0.038" J/Flex Guidewire
- 1-18G x 2 3/4" Introducer Needle.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733