

C₂ CryoBalloon®

180° Catheter

NEXT-LEVEL EFFICIENCY¹

**C2 CryoBalloon 180° Catheter for Cryoablation
of Long-Segment Barrett's Esophagus**

Delivers precise, targeted treatment while effectively
covering large areas.



C₂ CryoBalloon[®]

180° Catheter



An Innovation from C2 CryoBalloon—the 180° Catheter


Effectively treats long-segment Barrett’s esophagus while maintaining the flexibility, accuracy, and visibility of our focal catheter.

PRODUCT COMPARISON

	C2 CRYOBALLOON 180°	C2 CRYOBALLOON FOCAL	HALO 360°
Compliant balloon material	✓	✓	✗
Multiple balloon shapes	✓	✓	✗
Max. treatment length	3.0 cm	1.5 cm	4.0 cm
Min. treatment length	0.5 cm	1.5 cm	4.0 cm
Designed to specifically treat gastroesophageal junction adenocarcinoma (GEJ)	✓	✓	✗
Compatible with therapeutic scope	✓	✓	✗



180° CATHETER



FOCAL STANDARD CATHETER




A treatment for long-segment Barrett’s esophagus with proven reduction in post-treatment patient pain.²

- Lower post-treatment mean peak-pain score and duration of pain when compared to radiofrequency ablation (RFA)²

ORDERING INFORMATION

Catalog Number	Description
FG-1055	Standard 180° catheter
FG-1056	Pear 180° catheter

1. Compared to cryoablation using focal catheter. Data on file.
 2. van Munster et al. 2018. "Focal Cryoballoon versus Radiofrequency Ablation of Dysplastic Barrett’s Esophagus: Impact on Treatment Response and Postprocedural Pain." *Gastrointest Endosc* 88, no. 5 (Nov): 795–803.e2. doi: 10.1016/j.gie.2018.06.015.



Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
 1600 West Merit Parkway
 South Jordan, Utah 84095
 1.801.253.1600
 1.800.35.MERIT

Merit Medical Europe, Middle East & Africa (EMEA)
 Amerikalaan 42, 6199 AE
 Maastricht-Airport
 The Netherlands
 +31 43 358 82 22

Merit Medical Ireland Ltd.
 Parkmore Business Park West
 Galway, Ireland
 +353 (0) 91 703 733