


ReSolve® Thoracostomy Tray



ASSEMBLED TO INCREASE
PROCEDURAL **EFFICIENCY**

ASSEMBLED TO INCREASE PROCEDURAL EFFICIENCY.

The ReSolve Thoracostomy Tray has been designed with procedural ease in mind. Each tray component is placed in procedural order, helping to provide you with a smoother experience. To further support your procedures, a ReSolve catheter is included that offers several features designed to aid in correct placement.



RESOLVE[®] THORACOSTOMY TRAY INCLUDES:

Component	Quantity
+ Pigtail Drainage Catheter, 14F x 25 cm	1
+ Flexible Stiffening Cannula	1
+ ChloraPrep [™] Applicator (Orange Tint), 3 mL	2
+ Drape	1
+ Gauze, 4" x 4"	10
+ Lidocaine, 20 mL	1
+ Needle, 18G x 1.5"	1
+ Needle, 22G x 1.5"	1
+ Needle, 25G x 1.5"	1
+ Luer Lock Syringe, 10 mL	1
+ Luer Lock Syringe, 60 mL	1
+ Slip Syringe, 5 mL	1
+ Male-Female Non-vented Cap	1
+ Needle Access, 18G x 2"	1
+ Needle Access, 18G x 3.5"	1
+ Futura [®] Safety Scalpel, #11	1
+ ShortStop [®] Temporary Sharps Holder	1
+ Guide Wire, 0.038" x 80 cm	1
+ Dilator, 14F	1
+ 3-Way Stopcock	1
+ Connecting Tube	1
+ Heimlich Valve	1
+ Suture 0 Silk	1
+ Split Gauze, 4" x 4"	2
+ Clear Dressing	1

Large Drainage Holes


Achieve fast drainage with multiple large drainage holes.

Radiopaque Marker

Identify location of the most proximal drainage hole.

Depth Markers

Ideal reference for correct positioning and accurate placement.



ORDERING INFORMATION

Catalog Number	Product Description
RTT14038MB	ReSolve [®] Thoracostomy Tray

Packaged 4 units per box.

All trademarks are the property of their respective owners.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733