





EvenMore®

Chronic Hemodialysis Catheter

EvenMore Catheter Key Features


- Resistant Durathane® material
- Suture wing provides secure external anchoring
- Polyester cuff anchors catheter
- Ergonomic fittings
- Easy to read priming volumes printed on clamps



> Less than 5% recirculation, due to strategically placed venous return ports



> Luminal design allows for high flow



ORDERING INFORMATION

Catalog Number	Catheter Length	Description
EVENMORE 14.5F BASIC SET		
H787103035115†	20 cm	Straight Catheter
H787103035125	22 cm	Straight Catheter
H787103035015	24 cm	Straight Catheter
H787103035025	28 cm	Straight Catheter
H787103035035	32 cm	Straight Catheter
H787103035045	36 cm	Straight Catheter
H787103035055	40 cm	Straight Catheter
H787103035065	55 cm	Straight Catheter
H787103035075	24 cm	Pre-Curve Catheter
H787103035085	28 cm	Pre-Curve Catheter
H787103035095	32 cm	Pre-Curve Catheter
H787103035105†	36 cm	Pre-Curve Catheter

EVENMORE 14.5F CATHETER ONLY

H787103036015*	24 cm	Straight Catheter
H787103036025*	28 cm	Straight Catheter
H787103036035*	32 cm	Straight Catheter
H787103036065*	55 cm	Straight Catheter

†This item is made to order | Packaged 1 unit per box.

1 unit per carton. | *5 units per carton.

KIT CONFIGURATIONS INCLUDE:

BASIC SET

- | | |
|---|-----------------------------------|
| 1–EvenMore Dialysis Catheter | 2–Injection Caps |
| 1–15F Single Valve Peelable Introducer Sheath | 1–Adhesive Wound Dressing |
| 1–12F Vessel Dilator | 1–#11 Safety Scalpel |
| 1–14F Vessel Dilator | 1–0.038" J/Flex Guidewire |
| 1–Tri-Ball Tunneler with Sleeve | 1–18G x 2 3/4" Introducer Needle. |

Catalog Number	Description
H787065076025	15F EmboSafe™ Sheath Single Valve
H787055009012*	Marked Guidewire 3 mm J fixed core 0.038" x 70 cm uncoated depth markings every 5 cm
H787055009022*	Marked Guidewire 3 mm J fixed core 0.035" x 70 cm uncoated depth markings every 5 cm
H787103018015	Hemodialysis Catheter Repair Kit

5 units per carton. | *10 unit per carton.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733