




THREE LUMENS.  
ONE CATHETER.  
**PERFECT HARMONY.**

# trio-ct®

Triple Lumen Acute Hemodialysis Catheter

## INTERNAL LUMEN



17G lumen is independent from dialysis lumens and exits distally to them.

## CATHETER TIP



Tapered catheter tip design.

## KINK RESISTANT




13.5F oval catheter designed for kink resistance.

**Average recirculation rate of 2% in both the forward and reverse direction**

**Indicated for short term central venous access: hemodialysis, apheresis, infusion, power injection (PI) of contrast media, and central venous pressure (CVP) monitoring**

Up to 5mL/sec for PI of contrast media or a max pressure for PI of 300PSI



Rotating suture wings for secure external anchoring

Priming volume printed on catheter clamps

CVP monitoring and PI through third lumen

# THREE LUMENS. ONE CATHETER. **PERFECT HARMONY.**


## **POWER INJECTABLE LUMEN**

**Third lumen** decreases the need for additional sticks.

**Provides access for infusion** of medication, blood draws, central venous pressure (CVP) monitoring, and power injection (PI) of contrast media.

**Allows injection of contrast media** at a maximum of 5mL/sec.

**Extension is labeled** "Power Injectable."



## ORDERING INFORMATION

Catalog Number	Catheter Size	Catheter Length
<b>CATHETER ONLY—STRAIGHT</b>		
<b>ANG01135120</b>	13.5F	12 cm
<b>ANG01135150</b>	13.5F	15 cm
<b>ANG01135200</b>	13.5F	20 cm
<b>ANG01135240</b>	13.5F	24 cm
<b>ANG01135300</b>	13.5F	30 cm

Packaged 10 units per box.

Catalog Number	Catheter Size	Catheter Length
<b>BASIC KIT—STRAIGHT</b>		
<b>ANG01135121</b>	13.5F	12 cm
<b>ANG01135151</b>	13.5F	15 cm
<b>ANG01135201</b>	13.5F	20 cm
<b>ANG01135241</b>	13.5F	24 cm
<b>ANG01135301</b>	13.5F	30 cm

Packaged 5 units per box.

Catalog Number	Catheter Size	Catheter Length
<b>MAXIMAL BARRIER KIT—STRAIGHT</b>		
<b>ANG01135123</b>	13.5F	12 cm
<b>ANG01135153</b>	13.5F	15 cm
<b>ANG01135203</b>	13.5F	20 cm
<b>ANG01135243</b>	13.5F	24 cm
<b>ANG01135303</b>	13.5F	30 cm

Packaged 3 units per box.

## KIT CONFIGURATIONS INCLUDE:

### CATHETER ONLY KIT

- 1—Catheter
- 2—Dilators
- 1—Attachable Suture Wing
- 3—End Caps
- 1—Patient ID Card
- 1—Chart Sticker

### BASIC KIT

- 1—Catheter
- 2—Dilators
- 1—Attachable Suture Wing
- 3—End Caps
- 1—Patient ID Card
- 1—Chart Sticker
- 1—Introducer Needle
- 1—Guidewire
- 1—Scalpel,
- 1—Adhesive Wound Dressing

### MAXIMAL BARRIER KIT

- 1—Catheter
- 2—Dilators
- 1—Attachable Suture Wing
- 3—End Caps
- 1—Patient ID Card
- 1—Chart Sticker
- 1—Introducer Needle
- 1—Guidewire
- 1—Scalpel
- 1—Adhesive Wound Dressing
- 6—4" x 4" Gauze
- 1—Chlorhexidine
- 1—Tinted Chlorhexidine
- 3—10 mL Syringe
- 1—Filter Straw
- 1—Lidocaine/Saline Sticker Sheet\
- 1—3-0 Silk Curve Suture
- 2—Safety Hypodermic Needle
- 1—Hemostat
- 1—Full Body Max Barrier Drape
- 1—XL Gown
- 1—Bouffant Cap
- 1—Surgical Gloves
- 1—Needle Stick Pad
- 1—Steri-Strips
- 1—Face Mask
- 1—Surgical Face Mask

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

[merit.com](http://merit.com)

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733