




SCHON XL[®]

Acute Hemodialysis Catheter

SCHON XL[®]

Acute Hemodialysis Catheter



ORDERING INFORMATION

Catalog Number	Description
CATHETER ONLY	
H787108017015	15 cm Acute Catheter
H787108017025	20 cm Acute Catheter
H787108017035	24 cm Acute Catheter

Packaged 10 units per box.

Catalog Number	Description
CATHETER BASIC SET	
H787108007045	12 cm Acute Catheter
H787108007015	15 cm Acute Catheter
H787108007025	20 cm Acute Catheter
H787108007035	24 cm Acute Catheter

Packaged 5 units per box.

Catalog Number	Description
CATHETER FULL SET	
H787108027015	15 cm Acute Catheter
H787108027025	20 cm Acute Catheter
H787108027035	24 cm Acute Catheter

Packaged 5 units per box.

KIT CONFIGURATIONS INCLUDE:

BASIC SET

1-catheter
1-18G x 2 3/4" (7 cm) introducer needle
1-0.035" x 70 cm J/Flex guide wire
1-#11 blade disposable scalpel with handle
1-2-0 suture w/curved needle
1-HEMO-CATH[®] clip
1-adhesive wound dressing
1-12F vessel dilator
1-14F vessel dilator
2-injection caps

FULL SET

1-catheter
1-pair surgical gloves
4-4"x4" (10 cm) gauze sponges
1-fenestrated drape with clear adhesive
2-10 mL disposable luer lock syringes
1-15F Vascu-Sheath[®] Introducer
2-5 mL disposable luer lock syringes
1-25g x 5/8" (1.6 cm) needle
1-22g x 1 1/2" (3.8 cm) needle
1-straight disposable steel hemostat
1-5 mL ampule Lidocaine HCL 1%
1-18G x 1 1/2" (3.8 cm) needle
1-0.035" x 70 cm J/Flex guide wire
1-#11 blade disposable scalpel with handle
1-adhesive wound dressing
2-injection caps
1-18G x 2 3/4" (7 cm) introducer needle
1-HEMO-CATH clip
1-12F vessel dilator
1-14F vessel dilator
1-2-0 silk suture w/curved needle
1-3-pack Povidone Iodine swab sticks

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733