



THE GLOBAL LEADER IN WIRE-FREE LOCALIZATION



**TAKING THE
PRECISION OF SCOUT®
TO A NEW DIMENSION**



INTRODUCING SCOUT MD™

Improving surgical precision through refined RADAR technology.

The most widely adopted, clinically referenced wire-free platform transforms surgical localization with SCOUT MD™. Four distinct Reflector shapes and signals support:

- Precise tumor margin mapping (circumferential localization)
- Localizing multiple lesions in the breast, multiple nodes
- Localizing lesions only visible on MRI imaging




Detection mode available via Touch Screen
or SCOUTSelect on Guide

ADVANCED MARGIN AND TUMOR LOCALIZATION

Unique RADAR signatures facilitate marking multiple margins or lesions for greater surgical precision

ELIMINATE REDUNDANT PROCEDURES


The ONLY localization system that can be placed at the time of biopsy now provides 4 uniquely shaped markers.



SCOUT MD Reflector:

- Four distinct Reflectors can be differentiated on imaging (shapes) and in the OR (signals)
- Long-term implant status
- Clinically insignificant MRI artifact
- Can be placed anytime under the care continuum including prior to NAC and at the time of biopsy

SCOUT MD features two surgical guides to suit surgeon preference: Standard (18mm diameter) or Access (12mm diameter). The SCOUTSelect switch on both guides enables the surgeon to easily toggle between modes - to view distance readout of multiple reflectors or to isolate reflectors individually.



SCOUT MD Delivery Systems

- Can be placed under all imaging modalities (US, Stereo, CT, Mammo)
- All shapes available in 5cm, 7.5cm, 10cm and Single-hand



We are on a Mission to Reduce the Burden Cancer Places on Patients and Their Loved Ones



This brochure and other information are for the practitioner's convenience and for general information purposes only. This information does not constitute medical or legal advice, nor is it meant to endorse or guarantee the suitability of any of the referenced products or methods for any specific patient or procedure. Before using any product, refer to the Instructions for Use (IFU) for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle
East & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733