

# Rapid Response Chronic Dialysis Catheter Removal Pack

Designed  
alongside  
Physicians to  
help you reduce  
time and gain  
efficiencies



**EVERYTHING  
YOU NEED—  
EVERY TIME**

# Standard & Specialty Procedure Packs

Specializing in solutions that deliver efficiencies, Merit's Chronic Dialysis Catheter Removal Pack allows you to be prepared and in control every step of the way.

We've worked with interventionalists from across the country to develop a simple solution that includes all necessary equipment for a dialysis catheter removal procedure. This specialty pack can be implemented in your practice to save time, money, simplify ordering, reduce waste across your facility, and most importantly, allow you to focus on patient care.

## PACK BENEFITS\*

- Convenience
- Waste reduction
- High-quality componentry
- Cost control
- Supply consolidation
- Standardization of patient care
- Time savings



## CHRONIC DIALYSIS CATHETER REMOVAL PACK INCLUDES:

| Component Description                | Quantity |
|--------------------------------------|----------|
| + Transparent dressing, 4" x 4 3/4"  | 1        |
| + Needle, 25G x 1 1/2", reg bev      | 1        |
| + ChloraPrep™, 10.5 mL (clear)       | 1        |
| + L/L syringe, 10 cc                 | 1        |
| + Utility drape with tape, 15" x 26" | 1        |
| + Tray, 9" x 5" x 1"                 | 1        |
| + Scalpel, #15                       | 1        |
| + Curved iris scissor, 4 1/2"        | 1        |
| + Curved mosquito forceps, 5"        | 1        |
| + Sponge gauze, 2" x 2" (12-ply)     | 10       |
| + Sponge gauze, 4" x 4" (12-ply)     | 10       |

## AVERAGE COST SAVINGS\*

per procedure from implementing a catheter removal pack:

**\$3.00** Employee time-cost savings  
+ **\$2.00** Material cost savings  
+ **1.5** Minutes time savings

## AVERAGE YEARLY SAVINGS\*

Hospital: . . . . . **\$5,000**  
+ Ambulatory surgical center: . . . . . **\$29,000**  
+ Office-based lab . . . . . **\$2,000**

## ORDERING INFORMATION

| Catalog Number | Product Description                    | Units per box |
|----------------|--|---------------|
| K12T-12466B    | Chronic dialysis catheter removal pack | 10            |

Data on File. | Chloraprep is a trademark of Becton, Dickinson and Company.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.

1600 West Merit Parkway • South Jordan, Utah 84095 • Main: +1.801.253.1600 • Customer Service: +1.801.208.4300

merit.com

©2025 Merit Medical Systems, Inc. All rights reserved.

409052001\_002 ID 012125