

# SPEED MEETS PRECISION

## Embolization Gelatin Foam

Merit Medical embolization gelatin foam is uniformly cut and loaded into a syringe or torpedo cartridge for faster preparation and delivery.



S P E E D • P R E C I S I O N • C O N S I S T E N C Y™

# PRODUCT INFORMATION

Ready-to-use products, just hydrate with 50:50 saline/contrast mixture:

## EmboCube™ Embolization Gelatin



### EmboCube™ Embolization Gelatin

- Uniformly cut gelatin foam cubes preloaded into Medallion® syringe

## Torpedo™ Gelatin Foam



### Torpedo™ Gelatin Foam

- Preshaped torpedo gelatin foam loaded into a cartridge with optional blunt stylet

## EmboCube Embolization Gelatin Foam

| CATALOG NUMBER | HYDRATED SIZE | WEIGHT | RECOMMENDED INJECTION SYRINGE | MINIMUM CATHETER ID | COMPATIBLE MERIT PRODUCT                                                             | COLOR CODE |
|----------------|---------------|--------|-------------------------------|---------------------|--------------------------------------------------------------------------------------|------------|
| EC2550         | 2.5 mm        | 50 mg  | 1 mL                          | 0.024" (0.61 mm)    | Merit Maestro® Microcatheter 2.8F - 2.8F<br>Merit Maestro® Microcatheter 2.9F - 2.9F | ●          |
| EC25100        | 2.5 mm        | 100 mg | 1 mL                          | 0.024" (0.61 mm)    |                                                                                      | ●          |
| EC5050         | 5.0 mm        | 50 mg  | 3 mL                          | 0.040" (1.02 mm)    |                                                                                      | ●          |
| EC50100        | 5.0 mm        | 100 mg | 3 mL                          | 0.040" (1.02 mm)    |                                                                                      | ●          |

Packaged 5/box.

## Torpedo Gelatin Foam

| CATALOG NUMBER | HYDRATED SIZE | LENGTH | RECOMMENDED INJECTION SYRINGE | MINIMUM CATHETER ID | COMPATIBLE MERIT PRODUCT                 | COLOR CODE |
|----------------|---------------|--------|-------------------------------|---------------------|------------------------------------------|------------|
| TOR2510        | 2.5 mm        | 10 mm  | 1 mL                          | 0.027" (0.68 mm)    | Merit Maestro® Microcatheter 2.9F - 2.9F | ●          |
| TOR2520        | 2.5 mm        | 20 mm  | 1 mL                          | 0.027" (0.68 mm)    |                                          | ●          |
| TOR5010        | 5.0 mm        | 10 mm  | 3 mL                          | 0.040" (1.02 mm)    |                                          | ●          |
| TOR5020        | 5.0 mm        | 20 mm  | 3 mL                          | 0.040" (1.02 mm)    |                                          | ●          |

Packaged 5/box.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
Main: +1.801.253.1600  
Customer Service: +1.801.208.4300

Merit Medical Europe, Middle  
East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

703936001\_005 ID 101922