

Merit Advance[®]

ANGIOGRAPHIC NEEDLES

INTENDED USE

The Merit Advance[®] Angiographic Needle is intended for percutaneous vascular access.

CAUTIONS

- Read instructions prior to use.
- Store in a cool, dry place.
- This device is intended for single use only.
- Do not reuse or resterilize.

CONTRAINDICATIONS

There are no contraindications associated with the Merit Advance[®] Angiographic Needle.

INSTRUCTIONS FOR USE

1. Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
2. Insert a Merit Advance[®] Angiographic Needle using standard technique.

STORAGE CONDITIONS

Do not store in areas of excessive heat or humidity.

TRANSPORTATION CONDITION

Do not expose to excessive heat or humidity during transit.

SHELF LIFE

3 years



Manufacturer:
Merit Medical Systems, Inc.,
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.

Authorized Representative:
Merit Medical Ireland Ltd, Parkmore Business Park West,
Galway, Ireland

穿刺针及配件: Needles and Accessory

产品描述

用于经皮血管穿刺。

注意事项

- 使用前阅读说明书
- 贮存在低温干燥处
- 仅一次性使用
- 勿复用或再灭菌

禁忌症

配合使用造影穿刺针不会出现禁忌症状。

使用说明

1. 确定穿刺点，采用适当消毒技术处理穿刺点并按需求局部麻醉。
2. 采用标准技术扎入造影穿刺针。

贮存条件

勿贮存在高温或高湿处

运输条件

运输过程中无暴露在高温或高湿处

灭菌有效期

3年

型号和规格说明

型号	规格
AD18T51W	18G
AD18T71W	18G
AD18T71WT	18G
AD18N91W	18G
AD19U41W	19G
AD19U51W	19G
AD19U71W	19G
AD20T41W	20G
AD20T71W	20G
AD21T21W	21G
AD21T41W	21G
AD21T71W	21G
AD18T51WS	18G
AD18T71WS	18G
AD19U71WS	19G

AD18T71WC	18G
AD19U41WC	19G
AD19U71WC	19G
AD21T41WC	21G
AD21T91W	21G
SL18T71W	18G
SL18T71WS	18G
SL21T71W	21G
SL21T71WS	21G

注册证编号：国械注进20143036106

产品技术要求编号：国械注进20143036106

产品名称：穿刺针及配件

型号、规格：

AD18T51W, AD18T71W, AD18T71WT, AD18N91W, AD19U41W, AD19U51W, AD19U71W, AD20T41W, AD20T71W, AD21T21W, AD21T41W, AD21T71W, AD18T51WS, AD18T71WS, AD19U71WS, AD18T71WC, AD19U41WC, AD19U71WC, AD21T41WC, AD21T91W, SL18T71W, SL18T71WS, SL21T71W, SL21T71WS

结构及组成：

该产品由针管(304不锈钢)，针座(聚碳酸酯)和护套(低密度聚乙烯)组成。环氧乙烷灭菌，一次性使用，3年有效期。

灭菌方式：环氧乙烷灭菌

适用范围：该产品用于为血管介入器械进行经皮血管穿刺。

中国境内代理人及售后服务单位名称：麦瑞通医疗器械(北京)有限公司

中国境内代理人及售后服务单位住所：北京市朝阳区东大桥路9号楼2单元801室内B01、B02及B03单元。

中国境内代理人及售后服务单位电话：010-85610788

中国境内代理人及售后服务单位传真：010-85616981

注册人及生产企业名称：美国麦瑞通医疗设备有限公司

MERIT MEDICAL SYSTEMS, INC.

注册人及生产企业住所：1600 West Merit Parkway South Jordan, Utah 84095.USA

生产地址：1600 West Merit Parkway South Jordan, Utah 84095.USA

注册人及生产企业联系方式：1-801-253-1600

生产日期和失效日期：见产品标签。

说明书编制或修订日期：2021年9月