
C2 CryoBalloon Catheters Specifications


Balloon Shapes

- Catheter: Single Use, EtO sterilized,
- Compliant balloon material: Polyurethane (Latex Free)
- Nitrous Oxide temp -85 C
- Balloon pressure (~3.5 psig nom, 4.5 psig max)
- Procedure Requires Therapeutic Endoscope (T1 or T2) with 3.7mm or > channel size

Standard



Pear



C2 CryoBalloon Catheters Specifications

Diffuser Spray

Focal Spray Catheters

FG-1028: Standard Focal

FG-1024: Pear Focal

Focal Catheter Dosimetry:

Barrett's Esophagus: 8 sec

GAVE: 14 seconds



Diffuser Spray

180 Spray Catheters

FG-1055: Standard 180 spray

FG-1056: Pear 180 spray

180 Catheter Dosimetry:

Barrett's Esophagus: 1.2 mm/sec

